

Product datasheet for **MP221113**

Eml4 Mouse qPCR Primer Pair (NM_001114361)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	78798
Forward Sequence:	CTCTTGACGTGTGCTCAGGACA
Reverse Sequence:	CTGAGTGCGTTCCAATGGCTAC
ACCN:	BC065178 , BC066099 , BC067011 , BC070427 , BC079619 , NM_001114361 , NM_001114361.1 , BC065178.1 , BC032258 , BC051656
UniProt ID:	Q3UMY5
Synonyms:	4930443C24Rik; AI644019
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)