

## Product datasheet for **MP221128**

### Gpr160 Mouse qPCR Primer Pair (NM\_001134385)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	71862
Forward Sequence:	TGCCTTGCTTCAGGCTATGAGG
Reverse Sequence:	GGTCCAAGAAAGCACACGATG
ACCN:	<u><a href="#">BC116649</a></u> , <u><a href="#">NM_001134385</a></u> , <u><a href="#">NM_001134385.1</a></u> , <u><a href="#">NM_001134385.2</a></u> , <u><a href="#">BY036724</a></u> , <u><a href="#">BY208743</a></u> , <u><a href="#">BY527936</a></u>
UniProt ID:	<u><a href="#">Q3U3F9</a></u>
Synonyms:	1700025D19Rik; GPCR150
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)