

Product datasheet for **MP221137**

Fam110a Mouse qPCR Primer Pair (NM_146127)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	73847
Forward Sequence:	TCCAAGTCGGACTTGAGTGAGC
Reverse Sequence:	TAGAGACACGATGTCCGAGCTG
ACCN:	BC018162 , BC025092 , BC034240 , BC094463 , NM_146127 , NM_146127.1 , NM_146127.2 , NM_146127.3 , BY119831
UniProt ID:	Q8R184
Synonyms:	1700008J10Rik; 5430432M24Rik; A1747080
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)