

Product datasheet for **MP221176**

Lrrc8d Mouse qPCR Primer Pair (NM_001122768)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	231549
Forward Sequence:	CGACGCTGTCTTTGACCTCAC
Reverse Sequence:	GCTAAAGGCAGTCTGTCCACC
ACCN:	<u>BC037717</u> , <u>NM_001122768</u> , <u>NM_001122768.1</u> , <u>BC014756</u> , <u>BF148057</u> , <u>BM934577</u>
UniProt ID:	<u>Q8BGR2</u>
Synonyms:	2810473G09Rik; 4930525N13Rik; A930019F03; Lrrc5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)