

## Product datasheet for **MP221193**

### Chpf2 Mouse qPCR Primer Pair (NM\_133913)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	100910
Forward Sequence:	GCTTCAGTGCTCTGGAGTTGGA
Reverse Sequence:	AGCGAGAGTGTGGAGTAAAGG
ACCN:	<u><a href="#">BC047275</a></u> , <u><a href="#">NM_133913</a></u> , <u><a href="#">NM_133913.1</a></u> , <u><a href="#">NM_133913.2</a></u> , <u><a href="#">BC019714</a></u> , <u><a href="#">BC036970</a></u>
Synonyms:	2010209O12Rik; AW060945; mir-671; mKIAA1402
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)