

Product datasheet for **MP221259**

Pkig Mouse qPCR Primer Pair (NM_001039390)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 18769 |
| Forward Sequence: | AATGCAGTCCCTGACATCCAGG |
| Reverse Sequence: | AGGTGTTTCGATCACTGCTCTC |
| ACCN: | <u>BC026550</u> , <u>NM_001039390</u> , <u>NM_001039390.1</u> , <u>NM_001039390.2</u> , <u>BG086437</u> , <u>BY240874</u> , <u>BY266193</u> , <u>BY283174</u> , <u>BY343638</u> , <u>BY617655</u> |
| UniProt ID: | <u>O70139</u> |
| Synonyms: | PKIgamma |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)