

Product datasheet for **MP221371**

Snrpn Mouse qPCR Primer Pair (NM_001082961)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 20646 |
| Forward Sequence: | TTGGCACCTTCAAGGCTTTTGAC |
| Reverse Sequence: | GTAGCAAGACCAGACCCAAAACC |
| ACCN: | BC019589 , BC024880 , NM_001082961 , NM_001082961.1 , NM_001082961.2 , BC069894 , BQ042870 , BU614958 , BU708754 , BY246809 , BY269658 |
| UniProt ID: | P63163 |
| Synonyms: | 2410045I01Rik; HCERN3; Peg; Peg4; Pwc; sm-D; SMN; snRNP-N |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)