

Product datasheet for **MP221383**

Rgs3 Mouse qPCR Primer Pair (NM_134257)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	50780
Forward Sequence:	TGTCGGAGAACATCGCCAAGCA
Reverse Sequence:	GGAGAGGGATTTTGCTGGGTG
ACCN:	<u>BC033449</u> , <u>NM_134257</u> , <u>NM_134257.1</u> , <u>NM_134257.2</u> , <u>NM_134257.3</u> , <u>BC033449.1</u> , <u>BB320995</u> , <u>BB653029</u> , <u>BB739116</u> , <u>BC131680</u> , <u>BE686796</u> , <u>BY704855</u>
UniProt ID:	<u>Q9DC04</u>
Synonyms:	4930506N09Rik; C2pa; C2PA-RGS3; PDZ-RGS3; RGS3S
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)