

Product datasheet for **MP221446**

Cd247 Mouse qPCR Primer Pair (NM_001113394)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12503
Forward Sequence:	GCAGAAGCCTACAGTGAGATCG
Reverse Sequence:	CCAAGGGTCCTTGTGAGTTCTG
ACCN:	<u>BC052824</u> , <u>NM_001113394</u> , <u>NM_001113394.1</u> , <u>NM_001113394.2</u> , <u>BF020526</u> , <u>BM240376</u>
UniProt ID:	<u>P24161</u>
Synonyms:	4930549J05Rik; A430104F18Rik; AW552088; Cd3; Cd3-eta; Cd3-zeta; Cd3h; Cd3z; Cd3zeta; T3z; Tcrk
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)