

Product datasheet for **MP221650**

Kcnc1 Mouse qPCR Primer Pair (NM_001112739)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16502
Forward Sequence:	CGCTGGCAACCGCGCATCTG
Reverse Sequence:	GACAGAAGGTTGTGATGGAGACC
ACCN:	<u>BC125470</u> , <u>BC132439</u> , <u>NM_001112739</u> , <u>NM_001112739.1</u> , <u>NM_001112739.2</u> , <u>BB750064</u> , <u>BC100546</u> , <u>BE949887</u>
UniProt ID:	<u>P15388</u>
Synonyms:	C230009H10Rik; Kcr2-1; KShIIIB; Kv3.1; KV4; NGK2; Shaw
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)