

Product datasheet for **MP221654**

Enpp2 Mouse qPCR Primer Pair (NM_001136077)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18606
Forward Sequence:	GGATTACAGCCACCAAGCAAGG
Reverse Sequence:	ATCAGGCTGCTCGGAGTAGAAG
ACCN:	BC003264 , BC058759 , NM_001136077 , NM_001136077.1 , NM_001136077.2 , NM_001136077.3 , BB645672 , BY123373
UniProt ID:	Q9R1E6
Synonyms:	AT; ATX; Auto; E-NPP 2; lysoPLD; N; Npps2; Pd; PD-; PD-lalpha; Pdn2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)