

Product datasheet for **MP221711**

Pik3r1 Mouse qPCR Primer Pair (NM_001077495)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18708
Forward Sequence:	CAAACCACCCAAGCCCACTACT
Reverse Sequence:	CCATCAGCAGTGTCTCGGAGTT
ACCN:	BC026146 , BC051106 , NM_001077495 , NM_001077495.1 , NM_001077495.2 , BB647657 , BC002168 , BM963795 , BY276314
UniProt ID:	P26450
Synonyms:	p50alpha; p55alpha; p85alpha; PI3K
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)