

Product datasheet for **MP221718**

Camk2a Mouse qPCR Primer Pair (NM_177407)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12322
Forward Sequence:	AGCCATCCTCACCCTATGCTG
Reverse Sequence:	GTGTCTTCGTCTCAATGGTGG
ACCN:	<u>BC031745</u> , <u>NM_177407</u> , <u>NM_177407.1</u> , <u>NM_177407.2</u> , <u>NM_177407.3</u> , <u>NM_177407.4</u> , <u>BB659027</u> , <u>BY713276</u>
UniProt ID:	<u>P11798</u>
Synonyms:	CaMKII; mKIAA0968; R74975
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)