

Product datasheet for **MP221833**

Adnp2 Mouse qPCR Primer Pair (NM_175028)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	240442
Forward Sequence:	CATCGTGGGGAGATGTTTCTCC
Reverse Sequence:	GAGGGAAGTGAGCACCTTTGTTG
ACCN:	<u>NM_175028</u> , <u>NM_175028.1</u> , <u>BC024969</u> , <u>BC044898</u> , <u>BC044904</u> , <u>BC047352</u> , <u>BC064672</u>
UniProt ID:	<u>Q8CHC8</u>
Synonyms:	8430420L05Rik; BC024969; Kiaa0863; mKIAA0863; Zfp508
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)