

Product datasheet for **MP221946**

Slc4a4 Mouse qPCR Primer Pair (NM_001136260)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54403
Forward Sequence:	CAACGCCATCACTTTTGGAGGC
Reverse Sequence:	GCTGCTCAAGATGGTAAGCGGT
ACCN:	NM_001136260 , NM_001136260.1 , BC020415 , BC026592 , BC031756 , BC049236 , BC129863 , BC129864 , BC148292 , BC167254 , BF460877
UniProt ID:	O88343
Synonyms:	AI835705; NBC; NBC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)