

## Product datasheet for **MP221966**

### Hnrnpa1 Mouse qPCR Primer Pair (NM\_001039129)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15382
Forward Sequence:	CGAAACAACCGACGAGAGTCTG
Reverse Sequence:	CATGGCAGCATCCACTTCTCC
ACCN:	<a href="#">BC080675</a> , <a href="#">BC083136</a> , <a href="#">BC089340</a> , <a href="#">BC092395</a> , <a href="#">NM_001039129</a> , <a href="#">NM_001039129.1</a> , <a href="#">NM_001039129.2</a> , <a href="#">NM_001039129.3</a> , <a href="#">NM_001039129.4</a> , <a href="#">BC089340.1</a> , <a href="#">BC052723</a> , <a href="#">BE949469</a> , <a href="#">BY067931</a>
UniProt ID:	<a href="#">P49312</a>
Synonyms:	Hdp; HDP-1; hnrnp-A1; hnRNP A1; Hnrpa1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)