

## Product datasheet for **MP222021**

### Dst Mouse qPCR Primer Pair (NM\_134448)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13518
Forward Sequence:	GGTAACCTGCTGGAGACAGATG
Reverse Sequence:	CTGAGCCATGATCTTTCCGTGC
ACCN:	<u><a href="#">NM_134448</a></u> , <u><a href="#">BC019791</a></u> , <u><a href="#">BC033468</a></u> , <u><a href="#">BC048929</a></u> , <u><a href="#">BC131678</a></u> , <u><a href="#">BU919808</a></u>
UniProt ID:	<u><a href="#">Q91ZU6</a></u>
Synonyms:	2310001O04Rik; A830042E19Rik; ah; athetoid; AW554249; BP230; Bpag; Bpag1; BPAG1-n; dt; Macf2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)