

Product datasheet for **MP222049**

Use1 Mouse qPCR Primer Pair (NM_001145780)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67023
Forward Sequence:	CGTTACACCAGTGAGATGCGGA
Reverse Sequence:	TCGGATGCTGACTGCCTCTCAT
ACCN:	BC025086 , BC038519 , BC075695 , NM_001145780 , NM_001145780.1 , BC025086.1 , BY097923 , BY286283
UniProt ID:	Q9CQ56
Synonyms:	2010315L10Rik; 5730403H22Rik; AV002165; D12; Ed2; p31; Q-snare
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)