

Product datasheet for **MP222104**

Hdac6 Mouse qPCR Primer Pair (NM_001130416)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15185
Forward Sequence:	TCGCTGTCTCATCCTACCTGCT
Reverse Sequence:	GTCAAAGTTGGCACCTTCACGG
ACCN:	<u>BC041105</u> , <u>NM_001130416</u> , <u>NM_001130416.1</u> , <u>BC041105.1</u> , <u>BC033428</u> , <u>BM250542</u> , <u>BY195054</u> , <u>BY234769</u>
UniProt ID:	<u>Q9Z2V5</u>
Synonyms:	Hd6; Hdac5; mHDA2; Sfc6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)