

Product datasheet for **MP222118**

Pml Mouse qPCR Primer Pair (NM_178087)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18854
Forward Sequence:	GTCTAAGACCCAACCTGTGGCT
Reverse Sequence:	CTTCATGGAGCCGACTGTCTGA
ACCN:	BC020990 , NM_178087 , NM_178087.1 , NM_178087.2 , NM_178087.3 , NM_178087.4 , BB659153 , BG228200 , BX524753
UniProt ID:	Q60953
Synonyms:	1200009E24Rik; AI661194; Trim19
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)