

Product datasheet for **MP222119**

DbnDD2 Mouse qPCR Primer Pair (NM_026797)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	52840
Forward Sequence:	TCTTGCGGATCAGGATGGAGCA
Reverse Sequence:	GGATGTGGTCCTGTCAGACGTG
ACCN:	BC019947 , BC030298 , NM_026797 , BG080686
UniProt ID:	Q9CRD4
Synonyms:	1110017A21Rik; 2900022J10Rik; AU041050; AW048677; D2Bwg0891e; Nkip; R74724
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)