

## Product datasheet for **MP222195**

### Cant1 Mouse qPCR Primer Pair (NM\_001025617)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 76025  |
| Forward Sequence: | GGTCATCTACCAGATCGAAGGC   |
| Reverse Sequence: | CAGGTGTTTCGTCCTCACAGCA   |
| ACCN:             | <a href="#">BC020003</a> , <a href="#">NM_001025617</a> , <a href="#">NM_001025617.1</a> , <a href="#">NM_001025617.2</a> , <a href="#">BC023387</a> , <a href="#">BC048758</a> , <a href="#">BF465199</a>   |
| UniProt ID:       | <a href="#">Q8VCF1</a>   |
| Synonyms:         | 5830420C20Rik; Apy; Apy1h; D11Bwg0554e; Entpd8; SCA; SCAN-1; Sh; Shapy   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)