

## Product datasheet for **MP222198**

### Lrrc39 Mouse qPCR Primer Pair (NM\_175413)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 109245   |
| Forward Sequence: | GCTGAAAGAGAAGGTGACCAGG   |
| Reverse Sequence: | GGTCCTGTGAAGTTGCCATTCC   |
| ACCN:             | <u><a href="#">BC119529</a></u> , <u><a href="#">BC119530</a></u> , <u><a href="#">NM_175413</a></u> , <u><a href="#">NM_175413.1</a></u> , <u><a href="#">NM_175413.2</a></u> , <u><a href="#">NM_175413.3</a></u> , <u><a href="#">BC052542</a></u> ,<br><u><a href="#">NM_175413.4</a></u>  |
| UniProt ID:       | <u><a href="#">Q8BGI7</a></u>  |
| Synonyms:         | 2010005E21Rik; 9430028I06Rik; Myomasp  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)