

Product datasheet for **MP222329**

Gpx3 Mouse qPCR Primer Pair (NM_001083929)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14778
Forward Sequence:	ATCCTGCCTTCTGTCCCTGCTC
Reverse Sequence:	TGGTGAGGGCTCCATACTCGTA
ACCN:	BC003339 , BC037027 , BC049235 , BC061950 , NM_001083929 , NM_001083929.1 , BC037134 , BG796554 , BY061328 , BY100863 , BY339672
UniProt ID:	P46412
Synonyms:	AA960521; EGPx; GPx; GSHPx-3; GSHPx-P
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)