

## Product datasheet for **MP222335**

### Sh2b1 Mouse qPCR Primer Pair (NM\_011363)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20399
Forward Sequence:	ACAGCCCTAGAGATGCCTGACA
Reverse Sequence:	GCATTCCTGGATGTCAGACACC
ACCN:	<a href="#">BC011422</a> , <a href="#">BC051978</a> , <a href="#">NM_011363</a> , <a href="#">NM_011363.1</a> , <a href="#">NM_011363.2</a> , <a href="#">NM_011363.3</a> , <a href="#">BC011422.1</a> , <a href="#">BY749559</a>
UniProt ID:	<a href="#">Q91ZM2</a>
Synonyms:	AI425885; C530001K22Rik; Irip; mKIAA1299; Psm; SH2-B; SH2-Bb; Sh2b1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)