

Product datasheet for **MP222378**

Cep162 Mouse qPCR Primer Pair (NM_199316)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	382090
Forward Sequence:	CCAAGAAAAGCCTTCGCACCATG
Reverse Sequence:	AGCCTGGATGAGCTGTTACGCT
ACCN:	NM_199316 , NM_199316.1 , NM_199316.2 , BC151104 , BB413728 , BC004073 , BC028799 , BC064713 , BC150998 , BC151006 , BF465741 , BF721651
UniProt ID:	Q6ZQ06
Synonyms:	4922501C03Rik; mKIAA1009
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)