

Product datasheet for **MP222380**

Nadk2 Mouse qPCR Primer Pair (NM_001085410)

Product data:

| | |
|--------------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 68646 |
| Forward Sequence: | CGATGCTTCTGGCTGCGAGTAA |
| Reverse Sequence: | CGTGTATCGCACAGGCAGACAT |
| ACCN: | <u>BC115776</u> , <u>NM_001085410</u> , <u>NM_001085410.1</u> , <u>NM_001085410.2</u> , <u>BC048411</u> , <u>BE692075</u> , <u>BF460794</u> , <u>BG143421</u> , <u>BY545448</u> |
| UniProt ID: | <u>Q8C5H8</u> |
| Synonyms: | MNADK; Nadkd1 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |

