

Product datasheet for **MP222404**

Pard3 Mouse qPCR Primer Pair (NM_033620)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	93742
Forward Sequence:	GAGACTCTACGGAGGTCCATGT
Reverse Sequence:	TCGGTCATCCAGTTCTGTCTCG
ACCN:	BC090616 , NM_033620 , NM_033620.1 , NM_033620.2 , BC023413 , BC030315 , BC037126 , BC141187 , BC145508
UniProt ID:	Q99NH2
Synonyms:	AA960621; AI256638; Asip; D8Ertd580e; Par-3; Par3; Pard-3; Pard3a; Phip
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)