

## Product datasheet for **MP222488**

### Ctnnbip1 Mouse qPCR Primer Pair (NM\_001141930)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67087
Forward Sequence:	GCAAGGACAGTGGAGAATCAGC
Reverse Sequence:	TCATACCGTGCCTGCTCACACA
ACCN:	<u><a href="#">BC038253</a></u> , <u><a href="#">NM_001141930</a></u> , <u><a href="#">NM_001141930.1</a></u> , <u><a href="#">BC018436</a></u> , <u><a href="#">BQ031304</a></u> , <u><a href="#">BQ043263</a></u>
UniProt ID:	<u><a href="#">Q9JIN6</a></u>
Synonyms:	1110008O09Rik; 2310001I19Rik; AW457332; Ctnnbip1; Icat
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)