

Product datasheet for **MP222524**

Lamp2 Mouse qPCR Primer Pair (NM_010685)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16784
Forward Sequence:	GAGCAGGTGCTTTCTGTGTCTAG
Reverse Sequence:	GCCTGAAAGACCAGCACCAACT
ACCN:	<u>NM_010685</u> , <u>NM_010685.1</u> , <u>NM_010685.2</u> , <u>NM_010685.3</u> , <u>NM_010685.4</u> , <u>BC138718</u> , <u>BC138723</u> , <u>BQ176100</u> , <u>BY283193</u>
UniProt ID:	<u>P17047</u>
Synonyms:	CD107b; Lamp-2; Lamp-2a; Lamp-2b; Lamp-2c; Lamp II; LGP-B; Mac3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)