

Product datasheet for **MP222572**

Gprc5c Mouse qPCR Primer Pair (NM_001110337)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70355
Forward Sequence:	GGCTATGAGACCATCCTGAAGG
Reverse Sequence:	CAGCTGTCCATTGTAGCCACTG
ACCN:	BC008228 , BC031439 , NM_001110337 , NM_001110337.1 , BB871704 , BC005757 , BI134748 , BY003041 , BY108149
UniProt ID:	Q8K3J9
Synonyms:	1110028I06Rik; 3200002M13Rik; Raig3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)