

Product datasheet for **MP222612**

Mfap3 Mouse qPCR Primer Pair (NM_180599)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	216760
Forward Sequence:	GGAGTGTCTTCTCACGGTGGAT
Reverse Sequence:	CACGGTCATCAAAGGCTACGCT
ACCN:	<u>BC006828</u> , <u>NM_180599</u> , <u>NM_180599.1</u> , <u>BC058567</u> , <u>BY203406</u>
UniProt ID:	<u>Q922T2</u>
Synonyms:	2610509F16Rik; 2700079M14Rik; AW536327
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

