

Product datasheet for **MP300071**

Mir146b Mouse qPCR Primer Pair (MI0004665)

Product data:

| | |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | qPCR Primer Pairs |
| Forward Sequence: | GAGAACTGAATTCCATAGGC |
| Reverse Sequence: | GAACATGTCTGCGTATCTC |
| ACCN: | MIMAT0003475 |
| Synonyms: | MI0004665; MIMAT0003475; Mir146b mmu-mir-146b; mmu-miR-146b |
| Components: | 1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | The primer mix is stable for one year from date of shipping. Store at -20°C. |



[View online »](#)