

Product datasheet for **MP300562**

MIR1906-2 Mouse qPCR Primer Pair (MI0015984)

Product data:

Product Type:	qPCR Primer Pairs
Locus ID:	100316809
Forward Sequence:	TGCAGCAGCCTGAGGCAGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0007872
Synonyms:	MI0015984; MIMAT0007872; MIR1906-2 mmu-mir-1906-2; mmu-miR-1906
Components:	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)