

## Product datasheet for **MP300590**

### **MIR3067 Mouse qPCR Primer Pair (MI0014029)**

#### **Product data:**

|                   |                                                                                                                                                                                                                                                                                                                                                   |
|-------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                                 |
| Locus ID:         | 100526532                                                                                                                                                                                                                                                                                                                                         |
| Forward Sequence: | AGTTCTCAGGCCCGCTGTG                                                                                                                                                                                                                                                                                                                               |
| Reverse Sequence: | GAACATGTCTGCGTATCTC                                                                                                                                                                                                                                                                                                                               |
| ACCN:             | MIMAT0014840                                                                                                                                                                                                                                                                                                                                      |
| Synonyms:         | MI0014029; MIMAT0014840; MIR3067 mmu-mir-3067; mmu-miR-3067-5p                                                                                                                                                                                                                                                                                    |
| Components:       | 1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)                                                                                                                                                                                                                                                   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec. |
| Storage:          | The primer mix is stable for one year from date of shipping. Store at -20 °C.                                                                                                                                                                                                                                                                     |



[View online »](#)