

Product datasheet for **MP300644**

MIR344d-3 Mouse qPCR Primer Pair (MI0004227)

Product data:

Product Type:	qPCR Primer Pairs
Locus ID:	100526543
Forward Sequence:	AGTCAGGCTAGTGGTTATACT
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0014807
Synonyms:	MI0004227; MIMAT0014807; MIR344d-3 mmu-mir-344d-3; mmu-miR-344d-3-5p
Components:	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	The primer mix is stable for one year from date of shipping. Store at -20 °C.



[View online »](#)