

Product datasheet for **MP300657**

MIR3473b Mouse qPCR Primer Pair (MI0016997)

Product data:

Product Type:	qPCR Primer Pairs
Locus ID:	100628609
Forward Sequence:	GGGCTGGAGAGATGGCTC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0020367
Synonyms:	MI0016997; MIMAT0020367; MIR3473b mmu-mir-3473b; mmu-miR-3473b
Components:	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	The primer mix is stable for one year from date of shipping. Store at -20 °C.



[View online »](#)