

## Product datasheet for **MP300689**

### **MIR465d Mouse qPCR Primer Pair (MI0025028)**

#### **Product data:**

|                   |                                                                                                                                                                                                                                                                                                                                              |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                            |
| Forward Sequence: | TATTTAGAAATGGTACTGATGTG                                                                                                                                                                                                                                                                                                                      |
| Reverse Sequence: | GAACATGTCTGCGTATCTC                                                                                                                                                                                                                                                                                                                          |
| ACCN:             | MIMAT0029880                                                                                                                                                                                                                                                                                                                                 |
| Synonyms:         | MI0025028; MIMAT0029880; MIR465d mmu-mir-465d; mmu-miR-465d-5p                                                                                                                                                                                                                                                                               |
| Components:       | 1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)                                                                                                                                                                                                                                              |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | The primer mix is stable for one year from date of shipping. Store at -20°C.                                                                                                                                                                                                                                                                 |



[View online »](#)