

Product datasheet for **MP300744**

MIR5135 Mouse qPCR Primer Pair (MI0018047)

Product data:

Product Type:	qPCR Primer Pairs
Locus ID:	100628604
Forward Sequence:	AGGTCTAGGTGGCAAGGGCGT
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0020646
Synonyms:	MI0018047; MIMAT0020646; MIR5135 mmu-mir-5135; mmu-miR-5135
Components:	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)