

Product datasheet for **MP300874**

MIR669d-2 Mouse qPCR Primer Pair (MI0014051)

Product data:

Product Type:	qPCR Primer Pairs
Locus ID:	100526550
Forward Sequence:	ACTTGTGTGTCATGTATAT
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0005833
Synonyms:	MI0014051; MIMAT0005833; MIR669d-2 mmu-mir-669d-2; mmu-miR-669d-5p
Components:	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)