

Product datasheet for **MP300978**

MIR6993 Mouse qPCR Primer Pair (MI0022841)

Product data:

Product Type:	qPCR Primer Pairs
Forward Sequence:	AGTGGGAGAAACGGGTGTCTG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0027888
Synonyms:	MI0022841; MIMAT0027888; MIR6993 mmu-mir-6993; mmu-miR-6993-5p
Components:	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)