

Product datasheet for **MR200038L3V**

Pcp4 (NM_008791) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pcp4 (NM_008791) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pcp4
Synonyms:	P16Rimb19; Pcp-4; Pep19
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008791
ORF Size:	189 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200038).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008791.2
RefSeq Size:	669 bp
RefSeq ORF:	189 bp
Locus ID:	18546
UniProt ID:	P63054
Cytogenetics:	16 56.97 cM



[View online »](#)

Gene Summary:

Functions as a modulator of calcium-binding by calmodulin. Thereby, regulates calmodulin activity and the different processes it controls. For instance, may play a role in neuronal differentiation through activation of calmodulin-dependent kinase signaling pathways.
[UniProtKB/Swiss-Prot Function]