

## Product datasheet for **MR200046A1V**

### Mouse Hilpda (NM\_023516) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Hilpda (NM\_023516) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Hilpda  
**Synonyms:** 2310016C08Rik; AW049004; Hig2  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200046 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

**ATGAAGTTCATGCTGAACCTCTATGTGCTGGGCATCATGTTGACCCTGCTTTCCATCTTTGTTAGAGTGA**  
**TGGAGTCTCTGGGAGGCTTACTGGAGAGCCCACTGCCAGGAGCTCCTGGATCACGAGGGGTCAGCTAGC**  
**CAACACACAGCTCCCAAGGGCTGCCAGACCATCCATCCCGAGGAGTGCAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200046 protein sequence  
Red=Cloning site Green=Tags(s)

MKFMLNLYVLGIMLTLLSIFVRVMESLGLLESPLPRSSWITRGQLANTQPPKGLPDHPSRGVQ

**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_023516  
**ORF Size:** 195 bp  
**Buffer:** PBS with 0.001% Pluronic F68



[View online »](#)

**Stability:** AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

**RefSeq:** [NM\\_023516.3](#), [NP\\_076005.1](#)

**RefSeq Size:** 985 bp

**RefSeq ORF:** 195 bp

**Locus ID:** 69573

**UniProt ID:** [Q9JLS0](#)

**Cytogenetics:** 6 A3.3

**MW:** 7.1 kDa