

Product datasheet for **MR200046L3V**

Hilpda (NM_023516) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hilpda (NM_023516) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hilpda
Synonyms:	2310016C08Rik; AW049004; Hig2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_023516
ORF Size:	195 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200046).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_023516.3 , NP_076005.1
RefSeq Size:	985 bp
RefSeq ORF:	195 bp
Locus ID:	69573
UniProt ID:	Q9JLS0
Cytogenetics:	6 A3.3
Gene Summary:	Increases intracellular lipid accumulation. Stimulates expression of cytokines including IL6, MIF and VEGFA. Enhances cell growth and proliferation (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)