

Product datasheet for **MR200049L4V**

Uqcr10 (NM_197979) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Uqcr10 (NM_197979) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Uqcr10
Synonyms:	1110020P15Rik; AA960494; Ucrc
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_197979
ORF Size:	195 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200049).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_197979.2
RefSeq Size:	444 bp
RefSeq ORF:	195 bp
Locus ID:	66152
UniProt ID:	Q8R111
Cytogenetics:	11 A1
Gene Summary:	This is a component of the ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex), which is part of the mitochondrial respiratory chain. This subunit interacts with cytochrome c1 (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)