

Product datasheet for **MR200057A1V**

Mouse Gng13 (NM_022422) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Gng13 (NM_022422) AAV Particle
Tag: Myc-DDK
Symbol: Gng13
Synonyms: 1500031D04Rik; AB030194; Ggamma13
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200057 representing NM_022422
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGGAGTGGGATGTGCCCCAGATGAAGAAGGAGGTAGAGAGCCTCAAGTACCAACTGGCCTTCAAGA
GGGAGATGTCGTCCAAGACCATCCCCGAGCTTCTCAAGTGGATTGAGGATGGAATCCCAAGGACCCCTT
CCTGAACCCAGACCTGATGAAGAACAACCCTTGGGTAGAGAAGGCCAAGTGCACCATCCTA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200057 representing NM_022422
 Red=Cloning site Green=Tags(s)

MEEWDVPQMKKEVESLKYQLAFKREMSSKIPELLKWIEDGIPKDPFLNPDLMKNNPWVEKAKCTIL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_022422
ORF Size: 201 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [NM_022422.5](#), [NP_071867.1](#)

RefSeq Size: 438 bp

RefSeq ORF: 204 bp

Locus ID: 64337

UniProt ID: [Q9JMF3](#)

Cytogenetics: 17 A3.3

MW: 8.4 kDa