

Product datasheet for **MR200065**

Gng7 (NM_001038655) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Gng7 (NM_001038655) Mouse Tagged ORF Clone
 Tag: Myc-DDK
 Symbol: Gng7
 Synonyms: A1840417
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 ORF Nucleotide Sequence: >MR200065 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCAGGTAATAACAACGTCGCCAGGCCGGAAGCTGGTGGAGCAGCTGCGCATTGAAGCTGGGATCG
 AACGCATCAAGGTCTCCAAGCCTCGTCAGACCTGATGGGCTACTGTGAGCAACATGCCCGCAACGATCC
 CTTGCTGGTTGGCGTGCCAGCCTCTGAGAATCCATTCAAAGACAAAAGCCTTGATAATTCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200065 protein sequence
 Red=Cloning site Green=Tags(s)

MSGTNNVAQARKLVEQLRIEAGIERIKVSKASSDLMGYCEQHARNDPLLVGVPASENPFKDKKPCIIIL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001038655

ORF Size: 207 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001038655.1](#), [NP_001033744.1](#)

RefSeq Size: 3425 bp

RefSeq ORF: 210 bp
 Locus ID: 14708
 Cytogenetics: 10 39.72 cM
 MW: 7.5 kDa

Product images:



Circular map for MR200065