

Product datasheet for **MR200075A1V**

Mouse Crabp1 (BC065787) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Crabp1 (BC065787) AAV Particle
Tag: Myc-DDK
Symbol: Crabp1
Synonyms: AI326249; Crabp-1; CRABP-I; Crabpl; Rbp-5
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200075 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCCAACTTCGCCGTACCTGGAAGATGCGCAGCAGCGAGAATTTGACGAGCTCCTCAAGGCGTTGG
GTGTGAACCCATGCTGAGGAAGGTGGCCGTGGCGGCTGCGTCTAAGCCGCACGTGGAGATCCGCCAAGA
CGGGGATCAGTTCTACATCAAGACATCCACTACTGTGCGCACCCAGGAGATCAACTTCAAGGTCGGAGAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200075 protein sequence
 Red=Cloning site Green=Tags(s)

MPNFAGTWKMRSENFDELLKALGVNAMLKRVAVAAASKPHVEIRQDGDQFYIKTSTTVRTEINFKVGE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: BC065787
ORF Size: 210 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [BC065787](#), [AAH65787](#)

RefSeq Size: 709 bp

RefSeq ORF: 212 bp

Locus ID: 12903

Cytogenetics: 9 29.76 cM

MW: 7.8 kDa