

Product datasheet for **MR200093L3V**

2810047C21Rik (BC071238) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	2810047C21Rik (BC071238) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	2810047C21Rik
Synonyms:	2810047C21Rik; 100042870; AA683836; Gm4082
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC071238
ORF Size:	222 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200093).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC071238
RefSeq Size:	2423 bp
RefSeq ORF:	224 bp
Locus ID:	72716
Cytogenetics:	7 A1



[View online »](#)