

Product datasheet for **MR200094L4V**

Gngt1 (NM_010314) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gngt1 (NM_010314) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gngt1
Synonyms:	G(y)1; Gng1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010314
ORF Size:	225 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200094).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010314.1
RefSeq Size:	963 bp
RefSeq ORF:	225 bp
Locus ID:	14699
UniProt ID:	Q61012
Cytogenetics:	6 A1



[View online »](#)

Gene Summary:

Guanine nucleotide-binding proteins (G proteins) are involved as a modulator or transducer in various transmembrane signaling systems. The beta and gamma chains are required for the GTPase activity, for replacement of GDP by GTP, and for G protein-effector interaction. [UniProtKB/Swiss-Prot Function]